

Introduction

Key Benefits of the truCOVER[®] WGS Workflow

- Enhancing operational efficiency in sample preparation: integrated steps with or without amplification streamline the process, saving time and labor.
- Reliable DNA shearing across sample types and 2 orders of magnitude in input amounts without protocol adjustments.
- Lower total cost: fewer steps and consistent performance reduce repeats and save reagents and sequencing costs.




Figure 1. Overview of the truCOVER Workflow.

Methods

DNA Input and Shearing:

- High-quality gDNA (NA12878) and FFPE-derived DNA (FF12878) obtained from Coriell Institute
- FFPE DNA extracted using Covaris truXTRAC[®] FFPE SMART Kit on Hamilton Sonication STAR[®]
- Input range: 1–50 ng for library preparation
- Shearing performed using Adaptive Focused Acoustics[®] (AFA[®]) on Covaris R230 Focused-ultrasonicator




Library Preparation Workflow:

- Detailed protocols, including Covaris AFA settings, library quantification procedures¹, and workflow specifics are described in the truCOVER WGS Library Prep Kit product manual², prior publications^{3,4} and in the associated Application Note⁵

Library Quantification and Sequencing:

- All libraries were quantified using Qubit Fluorometer & pooled equimolar based on the calculated mass.
- Pooled libraries sequenced on Illumina NovaSeq[™] X Plus (~30× coverage) for performance verification.
- Sequencing data was analyzed using the CLC Genomics Workbench (Qiagen) Lightspeed Fastq to Germline Variant Pipeline.

Results



Conclusions

- Unlock Archival Samples:** Reliably generate whole-genome sequencing data from challenging FFPE samples to enable large-scale retrospective studies.
- Superior Variant Detection:** Achieves high F1 scores for both SNPs and INDELs from FFPE and high-quality DNA, ensuring accurate variant calling.
- Uniform Genome Coverage:** Delivers consistent, even coverage across regions of varying GC content, reducing analytical blind spots and improving data reliability.
- Conserve Precious Samples:** Requires minimal DNA input, preserving valuable and often irreplaceable FFPE tissue blocks for future studies.
- Complete Workflow Solution:** Benefit from a fully validated, single-vendor workflow from FFPE DNA extraction with truXTRAC SMART Solutions to library preparation and amplification with truCOVER DNA Library Prep Kit

References

- Ambavaram M., *et al.* (2025) Covaris truCOVER Library Prep Kit: A streamlined, PCR-free WGS workflow for enhanced NGS library performance and data quality. Covaris Application Note. https://www.covaris.com/wp/wp-content/uploads/resources_pdf/M020185.pdf
- Covaris. truCOVER DNA Library Prep Kits User Manual (2025). https://www.covaris.com/wp/wp-content/uploads/resources_pdf/010673_User_Instructions_truCOVER_DNA.pdf
- Process V., *et al.* (2025). Optimization of DNA Fragmentation Techniques to Maximize Coverage Uniformity of Clinically Relevant Genes Using Whole Genome Sequencing. *Diagnostics*, 15(18), 2294. <https://www.mdpi.com/2075-4418/15/18/2294>
- Amirault K., *et al.* (2025). Fully automated extraction of high-quality total nucleic acids from FFPE specimens for comprehensive genomic profiling of solid tumors. *SLAS Technology*, Volume 31, 100252 [https://www.slas-technology.org/article/S2472-6303\(25\)00010-X/fulltext](https://www.slas-technology.org/article/S2472-6303(25)00010-X/fulltext)
- Process V., *et al.* (2025) High-Fidelity Whole-Genome Sequencing from Low-Input, Formalin-Fixed Paraffin-Embedded (FFPE) DNA using the truCOVER[®] WGS Library Prep with Amplification. https://www.covaris.com/wp/wp-content/uploads/resources_pdf/mkt_5554.pdf

