


Introduction

Solving FFPE Challenges in Comprehensive Genomic Profiling (CGP)






- FFPE fixation introduces degradation and chemical damage that compromise DNA and RNA integrity, leading to uneven coverage, variable insert sizes, and reduced sensitivity for variant and fusion detection.
- Conventional FFPE workflows are typically single-analyte and rely on chemically fragmented nucleic acids, increasing variability, sample consumption, and performance loss in challenging genomic regions.
- We present an integrated, automated DNA-RNA co-extraction and library preparation workflow enabling high-fidelity multiomic profiling from FFPE tissue with improved coverage uniformity, variant calling accuracy, and fusion detection sensitivity.

Materials and Methods

- Samples:** Externally sourced FFPE scrolls (10–15 µm thick), n=48 for DNA extraction and WGS, n=28 for RNA-Seq
- FFPE extraction:** Automated truXTRAC[®] FFPE SMART Total NA (Hamilton[®] Sonication STAR[™]) vs Supplier A DNA/RNA FFPE kit
- WGS library preparation:** truCOVER[®] WGS Library Prep Kit with Amplification vs Supplier B Kit (100 ng input, 8 amplification cycles)
- RNA library preparation:** truCOVER Total RNA Kit vs Supplier C and Supplier D kits (50 ng input)
- Analysis:** Illumina[®] NovaSeq[™] X Plus (30× coverage, WGS); NextSeq[™] 2000 (20–25 million reads/sample, RNA-Seq); analysis via CLC Genomics Workbench and hap.py (DNA) and nf-core pipeline (RNA)



Results



Conclusions

- Integrated FFPE multiomics:** Automated DNA and RNA co-extraction (Hamilton Sonication STAR) enables parallel analysis from a single specimen reducing tissue consumption and variability.
- Superior DNA coverage performance:** Improved coverage uniformity supports reliable detection of SNPs and large INDELs (>16 bp) from degraded FFPE DNA.
- Variant calling accuracy:** Statistically significant improvements in precision, recall, and F1 score are observed for clinically actionable variants.
- RNA library quality:** Mechanical fragmentation yields consistent insert sizes, with residual rRNA (<1%) supporting improved gene expression analysis (higher exonic read content).
- Fusion sensitivity:** Complete detection of engineered gene fusions and splice variants is achieved at matched sequencing depth.
- Translational utility:** This complete workflow enables robust, high-fidelity multiomic profiling for cancer research and molecular discovery applications.
- Efficient end-to-end workflow:** Automated co-extraction and library preparation can be completed in ~18–20 hours.

References

- Arreaza, G., *et al.* (2016). Pre-analytical considerations for successful next-generation sequencing (NGS): Challenges and opportunities for formalin-fixed and paraffin-embedded (FFPE) tumor tissue samples. *International Journal of Molecular Sciences*, 17(9), 1579.
- Ambavaram, M., *et al.* (2025). Covaris truCOVER Library Prep Kit: A streamlined, PCR-free WGS workflow for enhanced NGS library performance and data quality. Covaris Application Note.
- McCarthy, P., *et al.* (2026). Whole transcriptome sequencing from degraded FFPE samples: A fast workflow to deliver robust gene fusion detection. Covaris Technical Note.
- Jaksik, R., *et al.* (2021). RNA-seq library preparation for comprehensive transcriptome analysis in cancer cells: The impact of insert size. *Genomics*, 113(6), 4149–4162.
- Amirault, K., *et al.* (2025). Fully automated extraction of high-quality total nucleic acids from FFPE specimens for comprehensive genomic profiling of solid tumors. *SLAS Technology*, 31, 100252.
- Process, V., *et al.* (2025). Optimization of DNA fragmentation techniques to maximize coverage uniformity of clinically relevant genes using whole genome sequencing. *Diagnostics*, 15(18), 2294.