


# Preparing High-Quality Proteins from FFPE Samples Just Got Easier



Designed specifically for high-quality protein sample preparation from formalin-fixed, paraffin-embedded (FFPE) samples, the Covaris truPREP® Protein FFPE Kit offers three modular options to facilitate easy implementation. This innovative solution provides the flexibility to employ either a comprehensive Covaris workflow or seamlessly integrate with your existing purification and/or digestion workflow. Streamline your sample preparation with easy loading of a 10 µm scroll into individually capped tubes for easy processing.

Achieve robust, reliable downstream results  
with the **truPREP Protein FFPE Kit**

Utilizing the power of focused ultrasonication with Adaptive Focused Acoustics® (AFA®) Technology, this workflow ensures:

- ✓ Complete tissue rehydration and emulsification of paraffin from FFPE samples (with lysis buffer).
- ✓ Consistent peptide yields (CV <20%) with no detergent carryover contamination.
- ✓ No degradation from toxic organic solvents or harsh chemicals.

Deparaffinization,  
Tissue Homogenization,  
& Decrosslinking



Reduction/Alkylation,  
Purification



Digestion



# Modular Options and Components

## Option 1 (PN 520359) Isolation

A proprietary non-toxic deparaffinization method using AFA on the R230 Focused-ultrasonicator effectively homogenizes tissues and extracts proteins. Customers who choose this option can continue with their established purification workflow.



## Option 2 (PN 520363) Isolation & Purification

Deparaffinization and extraction are combined with an optimized Protein Aggregation Capture (PAC) workflow for efficient protein purification. Customers who choose this option can utilize their established digestion workflow.



## Option 3 (PN 520363 & 520364) Isolation, Purification, and Digestion

This option includes the all three parts of the Covaris workflow: isolation, purification, and digestion. Promega Trypsin Gold is included for an overnight digestion, resulting in a purified peptide hydrolysate for further downstream analysis.



| Consumable/Reagent                          | Option 1 (PN 520359) | Option 2 (PN 520363) | Option 3 (PN 520363 & 520364) |
|---------------------------------------------|----------------------|----------------------|-------------------------------|
| prepTUBE TPX with Bead* (qty. 18)           | ✓                    | ✓                    | ✓                             |
| FFPE Protein Rehydration Buffer (1 mL   RT) | ✓                    | ✓                    | ✓                             |
| Deparaffinization Solution (8 mL   RT)      | ✓                    | ✓                    | ✓                             |
| Tissue Lysis Buffer (2 mL   RT)             | ✓                    | ✓                    | ✓                             |
| Collection Tube (qty. 16   RT)              |                      | ✓                    | ✓                             |
| Reducing Agent** (5 mg   -20 °C)            |                      | ✓                    | ✓                             |
| Alkylation Reagent (qty. 2   RT)            |                      | ✓                    | ✓                             |
| Magnetic Bead Mixture (100 µL   2–8 °C)     |                      | ✓                    | ✓                             |
| WB5 Diluent (2.4 mL   RT)                   |                      | ✓                    | ✓                             |
| Ethanol -Wash Buffer Diluent (1.6 mL   RT)  |                      | ✓                    | ✓                             |
| Promega Trypsin Gold*** (1 Vial   -20 °C)   |                      |                      | ✓                             |

RT - Room Temperature. \*Two additional prepTUBE TPX with bead have been included to be utilized during the Heat Block calibration step. These are not to be run with the assay after the calibration. \*\*Item can be stored at room temperature (RT) prior to reconstitution. Store at -20°C after reconstitution. \*\*\*Item is shipped at ambient temperatures but should be stored at -20 °C upon receipt of the product. Once reconstituted, the item should be used at once or, if only processing 8 samples at a time, the remaining material can be stored reconstituted at -20 °C.



14 Gill Street  
Woburn, Massachusetts  
01801-1721 USA



Tel: +1 781.932.3959



info@covaris.com

Covaris, truPREP, and AFA are registered trademarks of Covaris, LLC. All other trademarks are the property of their respective owners. Information subject to change without notice.

US: +1 781.932.3959  
customerservice@covaris.com

EU: +44 (0)845 872 0100  
emeacustomerservice@covaris.com

APAC: +86 137 6427 6714  
APACcustomerservice@covaris.com

Applications  
applicationsupport@covaris.com

Service and Support  
techsupport@covaris.com  
www.covaris.com | 2024© Covaris, LLC  
M020181\_RevA\_Aug2024

