



# truChIP<sup>®</sup> Ultra-Low Chromatin Shearing Kit

Adaptive Focused Acoustics (AFA)-based chromatin shearing for ChIP-based applications

Products PN 520156 and PN 520158

# Contents

| INTENDED USE  |
|---|
| INTRODUCTION  |
| REVISION HISTORY  |
| PROCEDURE OVERVIEW  |
| SAMPLE INPUT REQUIREMENTS                                     |
| KIT CONTENTS  |
| STORAGE   |
| SUPPLIED BY USER  |
| PROTOCOL  |
| A. Cell Preparation and Crosslinking6                         |
| A.1 Suspension cells  |
| A.2 Adherent cells  |
| B. Chromatin Shearing10                                       |
| SUPPLEMENTAL MATERIAL   |
| Appendix A: AFA Focused-ultrasonicator Operating Conditions12 |
| Appendix B: Chromatin Shearing Efficiency Analysis Protocol13 |
| Appendix C: Additional Notes                                  |
| REFERENCES  |

# **INTENDED USE**

The truChIP Kit is intended for use in research applications (RUO). This product is not intended for the diagnosis, prevention, or treatment of disease.

# INTRODUCTION

The truChIP Ultra-Low Chromatin Shearing Kit is optimized for the efficient and reproducible shearing of chromatin from adherent and suspension mammalian cell lines using Covaris AFA Focused-ultrasonicators. Focused-ultrasonicators provide a non-contact and isothermal method of shearing chromatin without compromising the structural integrity of the target epitopes of interest. This kit can be used to prepare sample for ChIP-qPCR, ChIP-ChiP, and ChIP-seq applications.

truChIP is compatible for use with a variety of cell lines, however, additional optimization for your specific cell line may improve results. To start, Covaris recommends for users to perform a one-time fixation and shearing time course study to empirically determine the optimal treatment conditions.

| Part Number | Revision | Date  | Description of change                             |
|-------------|----------|-------|---|
| 010255      | С        | 04/17 | Update template and publish ME220 settings        |
| 010255      | D        | 7/17  | Remove specific content information for Buffer D3 |
| 010255      | E        | 7/17  | Correct procedure overview                        |
| 010255      | F        | 11/17 | Added step in chromatin shearing protocol         |
| 010255      | G        | 1/19  | Update LE220 chromatin shearing settings          |

# **REVISION HISTORY**

#### **PROCEDURE OVERVIEW**

Collect cells and resuspend in fixing buffer Crosslink Protein-DNA interactions with formaldehyde Add Covaris shearing buffer 

Lyse and shear chromatin with an AFA® Focused-ultrasonicator

### SAMPLE INPUT REQUIREMENTS

The truChIP Ultra-Low Chromatin Shearing Kit is designed to perform efficient chromatin shearing of 100,000 cells or less.

| Single Sample                     | Low Cell      |
|-----------------------------------|---------------|
| Input cell number                 | <100,000      |
| Number of samples sheared per kit | 50            |
| AFA tube                          | microTUBE-130 |
| Shearing volume                   | 130 μl        |

#### **KIT CONTENTS**

| Buffer A     | 4 mL               | 10X Fixing Buffer                |
|--------------|--------------------|----------------------------------|
| Buffer D3    | 6 mL               | 10X SDS Shearing Buffer          |
| Buffer E     | 6 mL               | 1X Quenching Buffer              |
| Buffer F     | 0.8 mL             | 100X Protease Inhibitor Cocktail |
| Formaldehyde | 2x1 mL ampules (PN | 16% methanol-free formaldehyde   |
|              | 520156 only)       |                                  |

**Note:** Certain mammalian cell lines may have more proteases (nucleases), therefore, end-users can substitute **Buffer F** with other commercially available protease inhibitor cocktails if required.

Safety Data Sheets: http://covaris.com/resources/safety-data-sheets

## STORAGE

The kit is shipped cold and should be stored at 2-8C. Prior to use, kit reagent Buffers D3 and E may have to be warmed to 55C to dissolve precipitate and then cooled to room temperature before use.

Note: Mix solutions well before use to ensure solutions are completely solubilized

## **SUPPLIED BY USER**

- Molecular Biology Grade Water Thermo Scientific (Cat. No. SH3053802), Mo Bio (Cat. No. 17012-200), or equivalent
- Phosphate Buffered Salt Solution (PBS) Mo Bio (Cat. No. 17330-500), Thermo Scientific (Cat. No. SH30256.FS), or equivalent
- RNase A (DNase free) Thermo Scientific (Cat. No. EN0531) or equivalent
- Proteinase K (RNase and DNase free) Thermo Scientific (Cat. No. 17916), NEB (Cat. No. P8102S), or equivalent
- Covaris Focused Ultrasonicator (M220, ME220, S220, E220 Evolution, E220, or LE220)
- Refrigerated centrifuge having 15,000 x g capability
- Rocker Nutator<sup>®</sup> or equivalent
- AFA tubes & holders/racks

#### Ultra-Low Cell Protocol (<100K) Consumables & Holders/Racks Required

|        |                  | M220<br>Holder | ME220<br>Holder |          |         |          |        |
|--------|------------------|----------------|-----------------|----------|---------|----------|--------|
| Part   |                  | &              | &               | S-Series | E220Evo | E220     | LE220  |
| Number | Description      | Insert         | Insert          | Holder   | Rack    | Rack     | Rack   |
|        |                  |                |                 |          |         |          |        |
|        | microTUBE AFA    | 500414         | 500514          |          |         |          |        |
|        | Fiber Pre-Slit   | &              | &               |          |         |          |        |
| 520045 | Snap-Cap         | 500489         | 500526          | 500114   | 500433  | 500111   | NA     |
|        |                  |                | 500514          |          |         |          |        |
|        | microTUBE AFA    |                | &               |          |         |          |        |
| 520052 | Fiber Crimp-Cap  | NA             | 500526          | 500114   | 500433  | 500282   | 500282 |
|        | microTUBE-130    | 500414         | 500522          |          |         |          |        |
|        | AFA Fiber Screw- | &              | &               |          |         |          |        |
| 520216 | Сар              | 500489         | 500534          | 500339   | NA      | NA       | NA     |
|        |                  |                | 500514          |          |         |          |        |
|        | 8 microTUBE      |                | &               |          |         |          |        |
| 520053 | Strip V1         | NA             | 500526          | NA       | 500430  | 500191   | 500191 |
|        | 8 microTUBE-130  |                | 500518          |          |         |          |        |
|        | AFA Fiber Strip  |                | &               |          |         |          |        |
| 520217 | V2               | NA             | 500526          | NA       | NA      | NA       | NA     |
|        | 96 microTUBE     |                |                 |          |         | No rack  |        |
| 520078 | Plate            | NA             | NA              | NA       | NA      | required | 500329 |

# A. Cell Preparation and Crosslinking

Follow the Cell Preparation and Crosslinking method (A.1 – Suspension and A.2 – Adherent) for your cell culture type.

**Note:** ChIP assays are sensitive to crosslinking and shearing conditions. Therefore, Covaris recommends users to include multiple fixation and shearing time points to empirically determine the optimal treatment conditions.

# A.1 Suspension cells

1. Prepare solutions for the appropriate number of samples being processed fresh before starting.

| 8                  |   |
|--------------------|---|
| Buffer             | Instructions  |
| 1X Cold PBS        | Final Volume: 2.0 mL per sample                                       |
|                    | - Store on ice  |
| 1X Fixing Buffer A | Final Volume: 0.5 mL per sample                                       |
|                    | - Mix 50 μl of Fixing Buffer A with 0.450 mL of molecular biology     |
|                    | grade water   |
| Fresh 5%           | Final Volume: 1 mL per 1 to 20 samples                                |
| Formaldehyde       | - Mix 312.5 μl of 16% Fresh Formaldehyde with 687.5 μl of             |
|                    | molecular biology grade water   |
| Quenching Buffer E | Place in a 55C water bath to dissolve crystals, then place at ambient |
|                    |   |

**Note:** The use of fresh methanol-free formaldehyde is required to achieve reproducible results. The methanol-free formaldehyde ampule is for one-time use only–storage for later use is not recommended

2. Collect cells by centrifugation at 200 x g for 5 minutes at room temperature. Remove media and wash cells once with 1X Cold PBS and collect cells again by centrifugation.

| Reagent           | Ultra-Low Cell          |
|-------------------|-------------------------|
| 1X Cold PBS       | 400 μl                  |
| Input cell number | 1x10 <sup>5</sup> Cells |
| Centrifuge Tube   | 1.5 mL                  |

**Note:** Some cells do not pellet well at 200 x g. If a "spongy" pellet is not visible, increase speed at 100 x g intervals until a pellet is visible.

3. Resuspend cells in room temperature 1X Fixing Buffer A.

| Reagent         | Ultra-Low Cell |  |
|-----------------|----------------|--|
| Fixing Buffer A | 400 µl         |  |

4. Crosslink cells by adding freshly prepared 5% formaldehyde solution to a final concentration of 0.25% and set timer.

| Reagent               | Ultra-Low Cell |  |
|-----------------------|----------------|--|
| Fresh 5% Formaldehyde | 20 µl          |  |

**Note:** The use of fresh methanol-free formaldehyde is required to achieve reproducible results.

5. Place cells on a rocker at room temperature for the recommended time.

**Note:** We recommend including two fixation times. Typically, **2.5 and 5 minutes** for stem and primary cells, and **5 and 10 minutes** for all other cell types.

6. Quench the crosslinking reaction by adding the appropriate volume of Quenching Buffer E to the fixed cells. Keep cells on a rocker at room temperature for an additional 5 minutes.

| Reagent            | Ultra-Low Cell |  |
|--------------------|----------------|--|
| Quenching Buffer E | 12 μΙ          |  |

- 7. Collect cells by centrifuging at 500 x g for 5 minutes at room temperature.
- 8. Aspirate the supernatant and wash twice with cold PBS.

| Reagent     | Ultra-Low Cell |  |
|-------------|----------------|--|
| 1X Cold PBS | 300 μl         |  |

- 9. Collect cells by centrifugation at 500 x g for 5 minutes, 4C.
- 10. Proceed to nuclei preparation and chromatin shearing steps.

**Note:** You may flash-freeze the fixed cells in liquid nitrogen at this point and store at -80C for short periods of time (*e.g.*, 2 to 3 days). Longer-term storage is not recommended.

# A.2 Adherent cells

- 1. Grow the proper amount of cells to conduct a single ChIP assay or the initial time course until they are 80 to 90% confluent.
- 2. Prepare solutions for the appropriate number of samples being processed fresh before starting.

| Buffer  | Instructions   |  |
|---|--|--|
| 1X cold PBS   | Final Volume: 14 mL  |  |
|   | - Store on ice   |  |
| 1X Fixing Buffer A  | Final Volume: 2 mL   |  |
|   | - Mix 200 μl of Fixing Buffer A with 1.8 mL of molecular biology grade |  |
|   | water  |  |
| Fresh 5%  | Final Volume: 1 mL   |  |
| Formaldehyde  | - Mix 312.5 μl of 16% Fresh Formaldehyde with 687.5 μl of molecular    |  |
|   | biology grade water  |  |
| Quenching Buffer E  | Place in a 55C water bath to dissolve crystals, then place at ambient  |  |
|   |  |  |
| Important Notes   |  |  |
| The use of fresh methanol-free formaldehyde is required to achieve reproducible results |  |  |

- The methanol-free formaldehyde ampule is for one-time use only–storage for later use is not recommended
- 3. Remove media and wash each plate one time with cold PBS.

| Reagent | Volume |
|---------|--------|
| 1X PBS  | 2 mL   |

4. Remove PBS and add room temperature 1X Fixing Buffer A to each dish.

| Reagent         | Volume |
|-----------------|--------|
| Fixing Buffer A | 2 mL   |

5. Crosslink cells by adding freshly prepared 5% formaldehyde solution to a final concentration of 0.25% and start timing the crosslinking reaction.

| Reagent               | Volume |  |  |
|-----------------------|--------|--|--|
| Fresh 5% Formaldehyde | 100 μl |  |  |

**Note:** We recommend including two fixation times. Typically, **2.5 and 5 minutes** for stem and primary cells, and **5 and 10 minutes** for all other cell types.

6. Place cells on a rocker at room temperature for the recommended time.

**Note:** The use of fresh methanol-free formaldehyde is required to achieve reproducible results.

7. Quench the crosslinking reaction by adding the appropriate volume of Quenching Buffer E to the fixed cells. Keep cells on a rocker at room temperature for an additional 5 minutes.

| Reagent            | Ultra-Low Cell |  |  |
|--------------------|----------------|--|--|
| Quenching Buffer E | 60 µl          |  |  |

- 8. Completely aspirate the solution from the plate.
- 9. Add cold PBS to each dish and scrape cells from the plate into the proper vessel.

| Reagent         | Ultra-Low Cell |  |  |
|-----------------|----------------|--|--|
| 1X Cold PBS     | 450 μl         |  |  |
| Centrifuge Tube | 1.5 mL tube    |  |  |

10. Wash the plate with an additional volume of cold PBS to collect any remaining cells.

| Reagent         | Ultra-Low Cell |  |  |
|-----------------|----------------|--|--|
| 1X Cold PBS     | 450 μl         |  |  |
| Centrifuge Tube | 1.5 mL tube    |  |  |

11. Collect cells by centrifuging at 500 x g for 5 minutes, 4C.

**Note:** Some cells do not pellet well at 200 x g. If a "spongy" pellet is not visible, increase speed at 100 x g intervals until a pellet is visible.

#### 12. Wash cells twice by resuspending in cold PBS and collecting by centrifugation at 500 x g, 4C.

| Reagent         | Ultra-Low Cell |  |  |  |
|-----------------|----------------|--|--|--|
| 1X Cold PBS     | 450 μl         |  |  |  |
| Centrifuge Tube | 1.5 mL tube    |  |  |  |

- 13. Completely and carefully aspirate the supernatant from the tube(s) and place tube(s) on ice.
- 14. Proceed to chromatin shearing steps.

**Note:** You may flash-freeze the fixed cells in liquid nitrogen at this point and store at -80C for short periods of time (*e.g.*, 2 to 3 days). Longer-term storage is not recommended.

# **B.** Chromatin Shearing

1. Using the table below, prepare a sufficient volume of 1× Shearing Buffer D3 using the Buffer D3 and Buffer F stocks. A 15% excess volume is recommended when preparing this buffer.

| Total Number of<br>Samples | Buffer D3 – 10X<br>SDS Shearing<br>Buffer | Buffer F – 100X Protease<br>Inhibitor Cocktail | Molecular biology<br>grade water |
|----------------------------|---|--|----------------------------------|
| 1                          | 15 µl                                     | 1.5 µl   | 133.5 µl                         |
| 6                          | 90 µl                                     | 9 µl   | 801 µl                           |
| 12                         | 180 µl                                    | 18 µl  | 1.6 mL                           |
| 24                         | 360 µl                                    | 36 µl  | 3.2 mL                           |
| Х                          | X 15 µl                                   | X 1.5 µl                                       | X 133.5 µl                       |

\* Calculations include 15% excess

 Resuspend nuclei pellet in the Shearing Buffer D3 and transfer to appropriate AFA Tube(s). If conducting a shearing time course experiment, aliquot 130 µl of ≤100,000 fixed cells into 6 microTUBEs.

| Reagent            | Ultra-Low Cell |  |  |
|--------------------|----------------|--|--|
| Shearing Buffer D3 | 130 µl         |  |  |
| AFA Tube           | microTUBE-130  |  |  |

3. Shear chromatin with an AFA Focused-ultrasonicator with appropriate rack or holder; settings are provided in **Appendix A.** For the shearing time course, use processing times of 2, 4, 6, 8, 10, and 12 minutes.

**Note:** Optimization of shearing time should be conducted whenever experimental parameters (*e.g.,* cell type, cell number, or sample volumes) are changed.

- 4. After shearing, transfer samples into a pre-chilled microcentrifuge tube and place on ice until all tubes are processed. If batch processing using a high-throughput ultrasonicator (*e.g.* E220), samples can be maintained in the instrument's water bath at 4C before and after processing.
- 5. After processing, add 2-3 volumes of your IP buffer, and centrifuge samples at 10,000 x g, 4C for 5 minutes to pellet insoluble material

6. Transfer the supernatant to a new pre-chilled microcentrifuge tube.

**Note:** To check the efficiency of your shearing, reserve 130  $\mu$ l of the sheared chromatin and see **Appendix B** for detailed instructions

Note: Sheared chromatin can be stored at 4C for up to 2 days

**Note:** Freezing sheared chromatin is not recommended. Freeze/thaw cycles reduce IP efficiency and reproducibility

**Note:** For subsequent immunoprecipitation, sheared chromatin can be diluted in the desired immunoprecipitation buffer. Alternatively, the composition of the shearing buffer can be adjusted appropriately for immunoprecipitation. **The 1× SDS Shearing Buffer D3 composition is: 10 mM Tris-HCl pH 7.6, 1 mM EDTA, 0.1% SDS** 

# SUPPLEMENTAL MATERIAL

| Ultra-Low Cell Chromatin Shearing Protocol |          |          |            |                |          |          |
|--|----------|----------|------------|----------------|----------|----------|
| Instrument                                 | M220     | ME220    | S220       | E220 Evolution | E220     | LE220    |
| Target Size (bp)                           | 200-700  | 200-700  | 200-700    | 200-700        | 200-700  | 200-700  |
| PIP  | 75       | 75       | 105        | 105            | 105      | 300      |
| Duty Factor (%)                            | 5        | 5        | 2          | 2              | 2        | 15       |
| СРВ  | 200      | 1000     | 200        | 200            | 200      | 200      |
| Treatment Time (minutes)                   | 2-12     | 2-12     | 2-12       | 2-12           | 2-12     | 2-20     |
| Setpoint Temperature (C) $_1$              | 7        | 7        | 6          | 6              | 6        | 6        |
| Min/Max Temperature (C)                    | 4/10     | 6/12     | 3/9        | 3/9            | 3/9      | 3/9      |
| Sample Input                               | <100,000 | <100,000 | <100,000   | <100,000       | <100,000 | <100,000 |
| AFA Intensifier Required 2                 | NA       | NA       | Integrated | Yes            | Yes      | NA       |
| Water Level (run) 3                        | Full     | 9        | 12         | 6              | 6        | 6        |
| Sample Volume (µI) 4                       | 130      | 130      | 130        | 130            | 130      | 130      |

# Appendix A: AFA Focused-ultrasonicator Operating Conditions

#### Important Notes

1. If using the S220, E220 Evolution, E220, or LE220, set the temperature on the external chiller 3C below the setpoint temperature for the run. The min/max is set in SonoLab

2. If intensifier is required, please ensure PN 500141 is used

3. Water level should always be 1mm below the neck of the microTUBE-130 cap

4. Always fill the microTUBE-130 with 130 µl of sample

# Appendix B: Chromatin Shearing Efficiency Analysis Protocol

- 1. Take a 130  $\mu$ l aliquot of the sheared sample and transfer to 0.6 mL microcentrifuge tube.
- 2. Add 5  $\mu l$  of RNase A (10 mg/mL) and incubate at 37C for 30 min.
- 3. Add 5  $\mu$ l of Proteinase K (10 mg/mL) and reverse crosslink by heating at 65C overnight in a PCR cycler with a heated lid.
- 4. Purify DNA using either a commercial column based kit (*e.g.,* Qiagen QIAquick PCR Purification Kit, Cat. No. 28104), or phenol-chloroform extraction and ethanol precipitation.
- 5. Elute from column, or resuspend pellet with 25 μl of elution buffer (10 mM Tris-HCl, pH 8.5).
- 6. 1μL of purified DNA can be analyzed on an Agilent 2100 BioAnalyzer 12K chip to provide a more accurate representation of the shearing size range and distribution.
- 7. Alternatively, an aliquot of the sample can be run on a High Sensitivity Agilent 2100 chip.

# Appendix C: Additional Notes

- 1. The treatment settings listed in this document are recommended guidelines. Actual results may vary depending on the cell type and mass.
- 2. The Covaris process uses high intensity focused ultrasonic (HIFU) energy and as such is influenced by objects in the acoustic path from the transducer surface to the fluid sample. For example, particles and bubbles in the water bath may scatter the acoustic energy from the sample. Replace the bath water on a daily basis and ensure that appropriate time has been allowed for degassing and water bath temperature to stabilize prior to use of the instrument.
- 3. Bubbles in the sample fluid in the tube may diminish the acoustic dose effectiveness. Be sure to fill the tubes slowly with the recommended volumes and avoid the use of additional detergents that may induce foaming.

# REFERENCES

- 1. Sachs et al., Bivalent Chromatin Marks Developmental Regulatory Genes in the Mouse Embryonic Germline In Vivo, Cell Reports (2013)
- 2. Lee T.I., Johnstone S.E., Young R.A., Chromatin immunoprecipitation and microarray-based analysis of protein location. Nature Protocols (2006) 1:729-748.
- 3. Ralph M Bernstein, Ph.D. and Frederick C. Mills, Ph.D., Laboratory of Immunology, Division of Therapeutic Proteins, CDER, FDA,NIH Campus, Bethesda, MD. We very much appreciate their contribution to the shearing buffer SDS concentration titration experiment, formaldehyde fixation reduction time, and initial evaluation of our protocols and reagents.
- 4. Park P.J. ChIP-seq: advantages and challenges of a maturing technology. Nature Reviews Genetics (2009) 10: 669-680
- 5. Stewart D., Tomita A., Shi Y.B., Wong J., Chromatin immunoprecipitation for studying transcriptional regulation in Xenopus oocytes and tadpoles. Methods Mol Biol (2006) 322:165-182.
- Haring M, Offerman S, Danker T, Horst I, Peterhansel C and Stam M; Chromatin immunoprecipitation: optimization, quantitative analysis and data normalization, Plant Methods (2007), 3:11
- Mukhopadhyay A, Deplancke B, Walhout AJM and Tissenbaum HA; Chromatin Immunoprecipitation (ChIP) coupled to detection by quantitative real-time PCR to study transcription factor binding to DNA in *Caenorhabditis elegans*. Nature Protoc. (2008), 3(4) 698-70
- 8. Das P.M, Ramachandran K., vanWert J, Signal R.; BioTechniques (2004), 37:961-969