

A Robust LC-MS Analysis of *C. albicans* - Differential Protein Expression using Adaptive Focused Acoustics® (AFA®) Ultrasonication, SP3, and TMT Labeling

Covaris®

Authors: Nicolas Autret¹, Mathieu Cyrille², Xavier Méniche², Yannick Charretier², Leovigildo-Rey Alaban², and Frédéric Becquet²
Affiliation: 1 - Covaris, LLC., Woburn, Massachusetts 01801 and 2 - Bioaster, FR

Introduction: *Candida albicans* is the most common human pathogenic yeast and a member of the human gut flora [1]. It can reside harmlessly on the skin, in the oral cavity, and/or in the urogenital and gastrointestinal tracts causing cutaneous and mucocutaneous candidiasis [2]. It is a polymorphic organism, being able to grow in yeast form, or in hyphal and pseudohyphal types, which form filamentous structures [1]. *C. albicans* can cause life-threatening systemic infections: the hyphal form (which can be induced in vitro using BSA [3]) has an important role in causing disease by invading epithelial cells and causing tissue damage [4]. This study aimed at analyzing the protein determinants differentially expressed between yeast and hyphal forms.




Figure 1. The three morphological forms of *Candida albicans*. Pictures from Simon Vautier. <https://www.immunology.org/public-information/bitesized-immunology/pathogens-and-disease/candida-albicans>

Yeast Protein Extraction Workflow from Cell Lysis to TMT Labeling




Figure 2. Yeast Protein Extraction Workflow from Cell Lysis to TMT Labeling. Six SC5314 *C. albicans* isolates were cultured in 10 mL of either YPD medium (named "Yeast") or YPD medium + 10% BSA (named "Hyphae"). Proteins were extracted with AFA. Debris was removed by centrifugation and supernatant was transferred in a tube. Beads were added to perform clean-up and digestion through SP3 [5]. Peptides were labeled via TMT and pooled. LC-MS analysis was performed using Thermo Q Exactive High Resolution Accurate Mass Spectrometer.

Extraction with Precision - Adaptive Focused Acoustics (AFA) in Pre-analytical Sample Prep




Figure 3. Characteristics of AFA technology: A. Focalization of energy is permitted by high frequency acoustics dispensing. B. The resulting sonicated zone depends on the wavelength, corresponding to 1 to 3 mm in the case of AFA. C. Less power is needed with AFA, resulting in a non-contact and isothermal processing of the samples which ensures preservation of proteins, unlike conventional sonication techniques.

Protein and Peptide Quantification

All raw files that used DDA LC-MS/MS were analyzed by MaxQuant v1.6.0.1 software.
a - Proteins are quantified across 6 channels with ≥ 2 unique peptides and CV < 30% for each conditions.
b - Significantly differentiated proteins are defined as those with Benjamini-Hochberg adjusted p-value < 0.01 and a log₂ ratio exceeding ± 1 .

Proteins ^a	Unique Peptides	Total Peptides	Sig. Dif. Proteins Up	Sig Dif. Proteins Down ^b
2,993	29,148	29,648	65	130

1D SDS PAGE Analysis




Figure 4. 1D SDS PAGE analysis. Thirty µg of protein from each sample was analyzed on a 12% acrylamide SDS-PAGE gel (Biorad) under reducing conditions followed by Coomassie Brilliant Blue staining.

- Lane 1: MW marker
- Lanes 2 to 4: Yeast samples
- Lanes 5 to 7: Hyphae samples

Several bands seemed differentially expressed between the two phenotypes. The main observation came from a 60 kDa band clearly over-expressed in the Hyphae samples.

Differential Expression Between Yeast and Hyphae Subtypes




Figure 6. Differential expression between Yeast and Hyphae subtypes. Significantly differentiated proteins are defined as those with a Benjamini Hochberg adjusted p-value < 0.01 and a log₂ ratio exceeding ± 1 . In Hyphae, proteins in red are up-regulated and proteins in blue are down-regulated. Sixty five proteins are upregulated and 130 are down regulated in the Hyphae type. Those results are in line with previous observations. For example, Agglutinin-like protein 3 (ALS3, Q59L12), Hyphally regulated cell wall protein 1 (Q5AL03), and Induced during Hyphae Development Protein 1 (Q5A818) are known to be induced during Hyphae formation. Of particular interest, ALS3 is a key protein in the adherence of *C. albicans* hyphae to epithelial cells [6]. On the other hand, Yeast-form wall Protein 1 (Q59Y31) undergoes a downregulation under filamentation [7].

Principal Component Analysis of the Two Different Phenotypes




Figure 5. Principal Component Analysis of the two different phenotypes. The PCA based on the first and the second component, colored by strains shows a discrimination between Yeast and Hyphae on the first component with a variance of 84.3%. It also highlights the reproducibility of extraction, along with a clear separation between the two phenotypes.

Conclusion: This workflow optimally combines the most advanced preparation techniques for proteomics, namely the precise energy and thermal control of Covaris AFA with bead-based clean up (SP3) and TMT labeling. It can be easily adapted for automated, higher-throughput processing as already demonstrated for mammalian cells and tissues in 8-well strips [8] and 96-well plates [9].

References

- Gow NAR, Yadav B. Microbe Profile: *Candida albicans*: a shape-changing, opportunistic pathogenic fungus of humans. *Microbiology*. 2017;163(8):1145-1147. DOI: 10.1099/mic.0.000499
- Sudbery PE. Growth of *Candida albicans* hyphae. *Nat Rev Microbiol*. 2011;9(10):737-748. DOI: 10.1038/nrmicro2636
- Buu LM, Chen YC. Impact of glucose levels on expression of hypha-associated secreted aspartyl proteinases in *Candida albicans*. *J Biomed Sci*. 2014;21(1):22. DOI: 10.1186/1423-0127-21-22
- Meri T, Blom AM, Hartmann A, Lenk D, Meri S, Zipfel PF. The hyphal and yeast forms of *Candida albicans* bind the complement regulator C4b-binding protein. *Infect Immun*. 2004;72(11):6633-6641. DOI: 10.1128/IAI.72.11.6633-6641.2004
- Hughes CS, Moggridge S, Müller T, Sorensen PH, Morin GB, Krijgsvelde J. Single-pot, solid-phase-enhanced sample preparation for proteomics experiments. *Nat Protoc*. 2019;14(1):68-85. DOI: 10.1038/s41596-018-0082-x
- Nikou SA, Kichik N, Brown R, et al. *Candida albicans* interactions with mucosal surfaces during health and disease. *Pathogens*. 2019;8(2):53. DOI: 10.3390/pathogens8020053
- Granger BL, Flenniken ML, Davis DA, Mitchell AP, Cutler JE. Yeast wall protein 1 of *Candida albicans*. *Microbiology*. 2005;151(5):1631-1644. DOI: 10.1099/mic.0.27663-0
- Schreyer S, Berndt N, Eckstein J, et al. Dietary-challenged mice with Alzheimer-like pathology show increased energy expenditure and reduced adipocyte hypertrophy and steatosis. *Aging*. 2021;13(8):10891-10919. DOI: 10.18632/aging.202978
- Müller T, Kalxdorf M, Longuespée R, Kazdal DN, Stenzinger A, Krijgsvelde J. Automated sample preparation with SP 3 for low-input clinical proteomics. *Mol Syst Biol*. 2020;16(1). DOI: 10.15252/msb.20199111